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INSTRUCTION FOR USE

MODE D'EMPLOI

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Hearbow

für / for / pour / per / para / para

contact forte

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INTENDED USE: Hearbow is a wearing option designed to fulfill the intended purpose of the bone conduction hearing aids contact forte (BHM) and ADHEAR (med-el). Hearbow serves as a retention and coupling mechanism for these bone conduction hearing aids and ensures the transmission of energy from the bone conduction hearing aid to the patient without altering the energy (vibrations).

Hearbow is a non-invasive, non-active product.

USER GROUP(S) Hearbow is suitable for people of any age with hearing loss.

INDICATIONS: The indications for use of the compatible bone conduction hearing devices apply.

CONTRAINDICATIONS: Do not use in case of irritated and injured skin or pressure marks.

RESIDUAL RISKS: Possible pressure points, irritation of the skin or pain in case of wrong size selection of the wearing option.

CLEANING: Use a dry, clean cloth to clean the Hearbow. Make sure that your hands are also clean and dry. Remove sweat and moisture thoroughly, especially at points that are difficult to reach.

DISPOSAL: The packaging and parts of the Hearbow are designed to be disposed of in household waste. Please observe local regulations.

DECLARATION OF CONFORMITY: The full text of the declaration of conformity is available on request from BHM-Tech or at www.bhm-tech.com

**WARNINGS:**

- Hearbow may only be used on intact skin.
- Use of the Hearbow is only permitted after adjustment and detailed instructions have been given by your hearing aid acoustician or doctor.
- Please pay attention to the different sizes to ensure a comfortable fit.
Incorrect use and choosing the wrong size can lead to skin irritation, pressure points and pain.
- For safety reasons (risk of swallowing and choking), children or people with mental impairments may only wear the Hearbow under supervision.
- Take your Hearbow off when sleeping, as there is a risk of injury!
- Hearbow (contact forte) contains a magnet. Magnets can damage items such as credit and debit cards, data storage devices, mechanical watches, medical implants, etc. Keep magnets away from all devices and objects that could be damaged by magnetic fields.

**IMPORTANT INFORMATION:**

- All serious incidents occurring in connection with the product must be reported to the manufacturer of the product and to the competent authority.
- Only original accessories from BHM may be used. The product may only be used with the compatible hearing systems specified in these instructions for use.
- No alterations or modifications to the Hearbow are permitted!
- Defective devices must not be used. If the Hearbow is defective, please contact your hearing aid acoustician or doctor.
- If used correctly, the anticipated service life of the Hearbow is 12 months.

EXPLANATION OF SYMBOLS

.....

CE marking

Medical device

Attention! Observe the warnings
in the Instructions for Use

Follow the Instructions for Use!

Manufacturer

Date of manufacture

Reference number

Batch number

Unique device identifier

Temperature during transport
and storage

Ambient humidity during transport
and storage

For reuse on a
single patient

Protect against moisture

Prescription-only (only USA)

Important information about handling and product
safety

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